

## Staining Plants for GUS Activity

- Place plants (or organs) to be tested in a small Petri plate, a multiwell plate or a microfuge tube
- Cover the tissue with the GUS-staining medium (try to use as little as possible **EXPENSIVE**)
- Infiltrate tissue 3 times for 1 minute, each time “gently” releasing the vacuum
- Incubate 4-5 hours up to ON at 37°C
- remove media and discolor tissue by washing it 3 times for 90 minutes with 70% Ethanol
- assess staining macroscopically or microscopically

<b>Medium:</b>	372.5 ml	ddH <sub>2</sub> O
	7.5 ml	0.5M EDTA, pH 8.0
	3.2 ml	Triton X100 12%
	0.12 g	Sodium DiHydrogenoPhosphate (NaH <sub>2</sub> Phosphate)
	0.142 g	Di Sodium HydrogenoPhosphate (Na <sub>2</sub> H Phosphate)
	0.25 g	Potassium hexacyanoferrat (Fe <sup>3+</sup> )
	0.321 g	Potassium hexacyanoferric (Fe <sup>2+</sup> )

**Substrate** X-Gluc (5-bromo-4-chloro-3-indolyl glucuronide) **1mg/ml medium**